


1



Insert bolts in pill-shaped slots.
Then add washers and nuts.




3



Align upper stage
with bottom.

Camera lens should
be flush with stage
as shown.



2

This is an open hardware
design, developed by
hundreds of people
working collaboratively.




CERN OHL 1.2

* Visit PublicLab.org/micro for full instructions

Add washers and two more wingnuts to each
bolt as shown. Then thread rubber bands here.



4




5

Clip light on
and aim down
at slide

Public Lab

You're ready to
explore!



6

Place slide here

Adjust nuts at slots to

raise/lower stage & focus



Community Microscope Kit