


1

Insert bolts in pill-shaped slots.
Then add washers and nuts.



3

Align upper stage
with bottom.

Camera lens should
be flush with stage
as shown.

2

Attach converted
camera here by
clipping it
through this hole *




CERN OHL 1.2

This is an open hardware
design, developed by
hundreds of people
working collaboratively.

* Visit PublicLab.org/micro for full instructions

Add washers and two more wingnuts to each
bolt as shown. Then thread rubber bands here.

4



5

Clip light on
and aim down
at slide

Public Lab

You're ready to
explore!

6

Place slide here

Adjust nuts at slots to

raise/lower stage & focus



Community Microscope Kit